

Product Information

VACUETTE® URI-PLUS

Short Instructions

24-Hour Sample Collection



It is recommended that a separate clean cup or beaker be used to collect urine while voiding and the contents emptied into the 24h urine container each time.

Mix the contents of the 24h urine container prior to transferring a sample into the urine tube.

Sample Transfer



Lift up the label on the lid of the 24h urine container and firmly push the urine tube into the transfer device opening.

Hold the tube in position until the tube has filled completely.

For safety reasons, replace the label over the opening after specimen transfer.

Production location:
ENFA Envases Farmacéuticos, S.A., C/Paralela,15. Paracuellos de Jarama,
28860 Madrid, Spain
Distributed by Greiner Bio-One GmbH, Austria



Greiner Bio-One GmbH
Bad Haller Str. 32
A-4550 Kremsmünster



VACUETTE® URI-PLUS
24h urine container with integrated transfer device
For in vitro diagnostic use



Intended use

VACUETTE® URI-PLUS is to be used for collection of urine over a 24-hour period and subsequent transfer of a urine sample into urine tubes. The device is to be used by appropriately trained healthcare professionals in accordance with the instructions.

Product description

VACUETTE® URI-PLUS is a non-sterile, graduated, amber plastic container with an integrated transfer device in the lid including a long straw for sample collection into vacuum collection tubes.

Handling

The laboratory should give information on appropriate storage of urine during the 24-hour period.

Collection

1. After the patient has emptied his/her bladder into the toilet for the first time in the morning, the 24-hour collection begins. The time and date is noted on the container label and on any documentation provided as instructed by the laboratory.
2. During the following 24 hours, every drop of urine from the patient should be collected in the container. It is recommended that a separate clean cup or beaker be used to collect urine while voiding and the contents emptied into the 24h urine container each time.
3. Exactly 24 hours after beginning, the collection ends with a final emptying of the bladder into a cup or beaker to be added to the 24h urine container. Again, time and date should be documented on the label and on any associated paperwork.
4. The volume of urine collected should be determined by checking the fill against the scale on the container.

Transfer

5. Mix the contents of the 24h urine container prior to transferring a sample into the urine tube.
6. Lift up the label on the lid of the 24h urine container and firmly push the urine tube into the transfer device opening. Hold the tube in position until the tube has filled completely.
7. For safety reasons, replace the label over the opening after specimen transfer.

Cautions and precautions

- Do not touch the inside of the transfer device.
- Keep the label over the integrated transfer device during 24-hour collection to avoid accidental needlestick injury.
- Keep the container upright at all times to prevent any leakage of the contents.

Storage prior to collection

Keep away from direct sunlight.

Disposal

The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed. The lid contains a needle and should therefore be disposed of in an appropriate sharps container.

Labelling on packaging



Item number



LOT number



Temperature limit



Do Not Reuse



In Vitro Diagnostic Device



Manufacturer

Production location:
ENFA Envases Farmacéuticos, S.A., C/Paralela,15. Paracuellos de Jarama, 28860 Madrid, Spain
Distributed by Greiner Bio-One GmbH, Austria



Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria