



Tourniquets

from Greiner Bio-One

Practical helpers,
wide range

CONTENTS

1. **Product information**
Product line, Areas of application, Product advantages
2. **Application**
Single-use disposable tourniquet and Super-T
3. **FAQs**
4. **Quality**



PRACTICAL HELPERS
WIDE RANGE

PRODUCT INFORMATION

PRODUCT LINE

VACUETTE® SUPER-T
DISPOSABLE TOURNIQUET

VACUETTE®
TOURNIQUET
FOR CHILDREN

VACUETTE®
DISPOSABLE
TOURNIQUET

VACUETTE®
TOURNIQUET
(Standard product)


all GBO tourniquets do not contain
natural rubber latex

PRODUCT LINE

VACUETTE® TOURNIQUET

Standard product

- Re-usable
- Does not contain natural rubber latex
- Material:
 - Casing: Polycarbonate (PC)
 - Activation button: Polyoxymethylene (POM)
 - Material: Polyester/ Lycra



PRODUCT LINE

VACUETTE® TOURNIQUET FOR CHILDREN

- No pinching of skin on thin sensitive arms
- Children friendly print to distract and reduce anxiety
- Length: approx. 50cm
- Does not contain natural rubber latex
- With clip closure



PRODUCT LINE

VACUETTE® DISPOSABLE TOURNIQUETS


- Easy to place on patient's arm
- Quick release
- Protects against transfer of bacteria
- Single-use
- Dusted in talc for a comfortable application on the patient's arm
- no pinching



PRODUCT LINE

SUPER-T DISPOSABLE TOURNIQUET

- Contributes to reduction of nosocomial infection, e.g. by MRSA
- Easy handling
- Tear-resistant
- Soft silicon material
- Does not contain natural rubber latex, non-DEHP
- Lightly powdered
- Practical packaging



AREAS OF APPLICATION

Tourniquets are used in routine blood collection:

- To increase **VENOUS PRESSURE** in the arm
- To **IMPROVE PROMINENCE** of veins
- As well as make them easier to **FEEL** and **FIND**

Wherever there could be a risk of cross-contamination, especially with multi-resistant bacteria, it is recommendable to use disposable tourniquets.



PRODUCT ADVANTAGES

- **WIDE PRODUCT RANGE**
for various needs
- **DOES NOT CONTAIN NATURAL RUBBER LATEX**
- **EASY** handling
- Comfortable application on patients
- Practical **PACKAGING**

SUPER-T

- **TEAR-RESISTANT**
- **PROTECTION** against transfer of bacteria and infections (e.g. MRSA)
- **NON-DEHP**



DISPOSABLE
TOURNIQUETS

APPLICATION



1

Place the
tourniquet above
the puncture site.

Place one finger
underneath the
tourniquet at the
overlapping
position.



2

Begin to form a loop using the longer part of the tourniquet, pushing it down underneath the finger, through the rest of the tourniquet.



3

Pull the loop through.



4

Check that the
tourniquet is
positioned well.



5

Finished!





SUPER-T APPLICATION

You can find the complete instructions
for use on our website www.gbo.com

1

Place the
tourniquet above
the puncture site.



2

Pull the pointed end through the slit.



3

Pull until the
required tension is
achieved.



4

The tourniquet is now in position.

Blood collection can be proceeded with.



5

As soon as blood flows into the first tube, the tourniquet can be released with one hand.



Tourniquets

Application and products

FAQs



WHAT PRESSURE IS RECOMMENDED FOR APPLICATION

The ideal tourniquet pressure is 40mmHg.

HOW LONG SHOULD THE TOURNIQUET BE APPLIED FOR

One minute at most.

If left on for too long:

- Risk of haemolysis
- Inaccurate lab results
(Pumping with the fist can also distort results)

CAN SUPER-T BE USED MORE THAN ONCE

SUPER-T is a single-use product and should be disposed of immediately after use.

WHICH PATIENT GROUPS IS THE SUPER-T SUITABLE FOR

It is suitable for all adult patients.

There is no children's version
(with the position of the slit, that sometimes not enough
pressure can be created on smaller patients)



QUALITY AS SOLID BASIS



QUALITY AS SOLID BASIS

We have implemented a carefully thought-out **multi-level quality management system**.

- Production locations:
Certified according to **ISO 9001** and **ISO 13485**
- Products: **CE**-labelling
and **FDA**-registration

LEGAL REQUIREMENTS

according to medical device regulations all necessary information, which is required for a safe application of the device needs to be provided by the manufacturer.

Therefore please refer to the current valid instructions for use!

Download from: www.gbo.com/preanalytics

Greiner Bio-One tourniquets

PRACTICAL HELPERS.

WIDE RANGE



YOUR POWER FOR HEALTH