

# IN SALIVA VERITAS



GBO Saliva  
Collection System  
from Greiner Bio-One

# GBO Saliva Collection System for the

COLLECTION

of

SALIVA SAMPLES

from both adults and young persons

for

ORAL FLUID

TESTING

# Oral Fluid Testing

DRUGS OF ABUSE

PERIODONTITIS

DNA – LIFESTYLE  
(Nutrigenetics e.g. Lactose intolerance)

ORAL CANCER  
(HPV, cancer metabolites, RNA)



# DRUGS OF ABUSE

- Clinical drug testing
- Workplace drug testing
- Legal medicine - prisons
- Roadside testing
- Antidoping agencies



# GBO Saliva Collection System

The first system on  
a liquid basis:

## GBO Saliva Collection System (SCS)





Collection within 3 minutes: 1-2-3

# INSTRUCTIONS

For the complete handling recommendation  
please refer to the Instructions for Use  
for each product on [www.gbo.com](http://www.gbo.com).



# Collection within 3 minutes: 1-2-3

For the complete handling recommendation please refer to the Instructions for Use for each product on [www.gbo.com](http://www.gbo.com).



1

## Oral Fluid Collection

Rinse oral cavity  
**2 MIN** with **SALIVA  
EXTRACTION  
SOLUTION N°1**

2

## Saliva Collection Beaker

Unscrew **SALIVA  
COLLECTION BEAKER  
N°2** and spit in entire  
contents of mouth.

3

## Transfer

Screw lid gently back onto beaker 2, remove safety sticker and take **SALIVA TRANSFER TUBE N°3** and push down onto the opening. Tube will fill by itself.

## STORAGE OF SAMPLES

Depending on the analytes

**+20 to +25°C**  
storage 1 day - 20 days

**+4°C**  
storage 1 week - 3 weeks

**-20°C**  
storage for months

## CENTRIFUGATION

Saliva Transfer Tubes  
containing saliva samples  
should be centrifuged

for **10 MIN** at **2200g**.





For the determination  
of the saliva content in  
the collected sample:

## SALIVA QUANTIFICATION KIT

### 5 CALIBRATORS

with defined saliva content:  
11, 27, 46, 66 and 85 U/mL

correspond to 11 – 85 volume  
percentage of saliva

### 2 CONTROLS:

30 U/mL and 70 U/mL




# IN SALIVA VERITAS

The yellow food dye tartrazine, which is part of calibrators, controls and collected samples, enables to determine the content **PHOTOMETRICALLY**.

Using the calibration line, the saliva content (U/ml) in the sample is calculated using the linear regression formula.

## Calibration and Saliva Quantity Determination

$$\text{Extinction cal} = \text{Ex } 450 \text{ nm} - \text{Ex } 520 \text{ nm}$$



# IN SALIVA VERITAS

Amount (%) of  
oral fluid important  
for re-calculation  
of the analytes:

*e.g. Cocaine, cut off 10 ng/mL*

COCAINE MEASURED	% OF ORAL FLUID IN SAMPLE	COCAINE IN 1 ML ORAL FLUID
5 ng	58 %	8,6 ng
		negative
5 ng	45 %	11,1 ng
		positive



Simple, safe, fast.

# BENEFITS



# Benefits

**NO TIME DELAY** caused by dry mouth  
(fixed sampling time of 2 minutes)

pH 4.2 of the saliva extraction solution  
**STABILIZES ANY PARENT DRUG  
AND METABOLITE** immediately

**NO SAMPLE ADULTERATION**  
by pressing or sucking on the carrier  
material – no loss of analytes

Internal standard (Tartrazine) enables  
**DONOR SPECIFIC QUANTIFICATION**

**WHOLE RINSING SAMPLE**  
can be collected



# Benefits

pH stays **ACIDIC** (pH ~5.5)  
6AM, Cocaine, Z-substances ...

No detergents and no ion-suppression  
on **MASS SPECTROSCOPY**  
INSTRUMENTS

**IMMEDIATE SAMPLE SPLITTING**  
in two or more sample containers

**SAMPLE ADULTERATION EASILY DETERMINED** by measurement of  
Amylase (EIA) and Cortisol (MS)

**FILL VOLUME IS READABLE**  
– no sucked pads have influence





# Benefits

**EASY PROCESS INTEGRATION**  
and short turn-around time

**STANDARDISED TUBES**  
for all racks and analyzers

**CLOSED SYSTEM** for hygiene  
and sample integrity

Analysis can be conducted  
**FROM THE ORIGINAL  
SAMPLE CONTAINER**

Clearly arranged **STORAGE,**  
**AUTOMATION** also possible



# Oral Fluid Testing

DRUGS OF ABUSE

PERIODONTITIS

DNA – LIFESTYLE

(Nutrigenetics e.g. Lactose intolerance)

ORAL CANCER

(HPV, cancer metabolites, RNA)

- **HUMAN GENES (SNP)**  
e.g. Lactose intolerance
- **ORAL CANCER**  
(HEAD & NECK CANCER)  
Human or cancer cell DNA, mRNA  
and cancer metabolites
- **PERIODONTITIS**  
periodontitis associated pathogenic  
bacterias: ParoCheck®
- **CLINICAL DIAGNOSTIC  
OF INFECTIONS**  
direct detection of viral pathogens  
e.g. HSV I/II or VZV-DNA

## MOLECULAR DIAGNOSTIC (PCR)

Detection of  
DNA und RNA  
with molecular  
biological methods

The system is made up  
of various components.

# PRODUCT OVERVIEW



# GBO Saliva Collection System

- For self collection (single packed) [Item N° 881001](#)
- For professional use [Item N° 881002](#)  
50 Sets with all components in bulk  
GBO Saliva Extraction Solution  
GBO Saliva Collection Beaker  
GBO Saliva Transfer Tubes
- For professional use  
GBO Saliva Extraction Solution [Item N° 881112](#)  
GBO Saliva Collection Beaker [Item N° 881210](#)  
GBO Saliva Transfer Tubes [Item N° 881113](#)  
(with sleeves)  
or GBO Saliva Transfer Tubes [Item N° 881117](#)  
(without sleeves)



# GBO Saliva Quantification Kit

- For professional use in laboratories  
Item N° 881010







# QUALITY AS SOLID BASIS



# QUALITY AS SOLID BASIS

We use a carefully designed multi-level quality assurance system.

- Production locations:  
ISO 9001 and ISO 13485  
certified
- Products are CE certified

# YOUR POWER FOR HEALTH