

Spinner Flasks with IOB Impeller, Standard and Water Jacketed

- Internal impeller decreases risk of contamination through cap
- Adjustable shaft aligns paddle with bottom of flask
- Large paddle promotes increased oxygenation with low shear



Stock No.	Description	Center Cap (mm)	Sidearm Cap (mm)	OD x OAH (mm)
1965-61001	Complete, 100mL	70	32	65 x 135
1965-61002	Complete, 250mL	70	32	85 x 140
1965-61005	Complete, 500mL	70	45	110 x 180
1965-61015	Complete, 500mL	100	45	110 x 180
1965-61010	Complete, 1L	100	45	135 x 250
1965-61030	Complete, 3L	100	45	165 x 280
1965-61060	Complete, 6L	100	45	204 x 410
1965-61080	Complete, 8L	100	45	237 x 430
1965-61150	Complete, 15L	100	45	294 x 460
1965-61360	Complete, 36L	100	45	410 x 535

Water Jacketed Spinner Flasks

1965-62001	Complete, 100mL	70	18	65 x 135
1965-62002	Complete, 250mL	70	24	85 x 140
1965-62005	Complete, 500mL	70	38	110 x 180
1965-62010	Complete, 1L	100	45	135 x 250
1965-62030	Complete, 3L	100	45	165 x 280
1965-62060	Complete, 6L	100	45	204 x 110

Internal Overhead Bearing Impeller with Cap

- Stainless steel shaft design prevents cell damage and minimizes magnetic decoupling
- Sealed internal assembly requires less incubator space
- For use with Bellco 1965-series Standard and Water Jacketed Spinner Flasks
- 3L and up include adjustable paddle blades



Impeller Shown with adjustable paddle blade.

Stock No.	Paddle L x W (in)	Flask (mL)	Cap (mm)
1965-60001	2 x 1.5	100 & 250	70
1965-60015	2 x 2	500	70
1965-60005	2 x 2	500	100
1965-60010	4 x 3	1000	100
1965-60030	4 x 3	3000*	100
1965-60060	6.75 x 5	6000*	100
1965-60080	6.75 x 5	8000*	100
1965-60150	9.5 x 6	15000*	100
1965-60360	9.5 x 6	36000*	100

*Includes adjustable width paddle