

CELLSTAR® Cell Culture Microplates

Cell culture treated microplates are available in the following versions: 96, 384, 1536 well format

Properties

- Improved cell adhesion through physical surface treatment
- Compatible with automated systems. For further information concerning ANSI standards please visit our website: www.gbo.com/bioscience/technical_information
- Alphanumeric well coding



An overview of all 96 well, 384 well and 1536 well microplates listed in this catalogue can be found in the Technical Appendix → p. A 1 3 ff.



Detailed technical information on microplates
→ p. 2 1 2 ff. HTS microplates
→ p. 3 1 2 ff. Immunology
→ p. 14 1 4 Barcode labelling of microplates

Further literature on cell culture treated microplates

- **Application Note “Establishing a cell culture assay based on TR-FRET for screening G-Protein-coupled receptors”** (F074 058)
- **Application Note “Selection of cell culture surfaces for the adipogenic differentiation of hMSCs”** (F010 003)

96 Well Polystyrene Cell Culture Microplates

Cell culture treated 96 well microplates are available in the following versions:

- With U-, V- and F-bottom
- Clear, black and white
- Standard or half area microplates
- Black and white “clear bottom” plates (μ Clear®)
- Chimney well design, raised wells and condensation rings in lids prevent cross-contamination
- With or without lid
- Improved cell adhesion through physical surface treatment
- Barcode-labelled on request

Properties

- Lid enables gas exchange with minimal evaporation
- High clarity of the clear microplates for optimal microscopic examinations
- Stackable
- Alphanumeric well coding
- Individually wrapped peel-off bags
- Consecutive lot numbering

Well Profile

Depending on the application, the well profile is a key feature in a 96 well cell culture microplate. For further information and figures on the well profiles see → p. 2 1 6 f.

1. U-Bottom

The “U” describes the round bottom shape.

- No sharp corners to facilitate easy and residue-free pipetting
- Suitable for +/- analyses
- Working volume: 40 – 280 μ l

2. V-Bottom

The “V” stands for the conically tapered well bottom.

- For precise pipetting
- Suitable for +/- analyses
- Working volume: 40 – 200 μ l

3. F-Bottom / Standard (ST)

The “F” refers to the flat well bottom.

- Excellent optical properties
- For precise optical measurements
- For microscopic applications (bottom reading)
- Cell growth area: 32 mm²
- Working volume: 25 – 340 μ l

4. F-Bottom / Chimney Well

The chimney well cell culture microplate has the same well profile as the standard F-bottom plate. The difference to the standard plate is the chimney-like arrangement of the wells i.e. each well stands on its own. Therefore the risk of contamination from sample material being carried over is minimised.

- Cell growth area: 34 mm²
- Working volume: 25 – 340 μ l

μ Clear® / Solid Bottom

Clear bottom microplates have pigmented walls and a transparent thin film bottom, the so-called μ Clear® bottom. In contrast to our standard microplates with a solid polystyrene bottom, they are ideal for cell culture and microscopic applications using fluorescence or luminescence detection methods.

Half Area Microplates

For many applications, a reduction of the sample volume is an important feature. Beside high-format plates, the 96 well half area microplates offer an interesting alternative here. They can be pipetted automatically as well as manually without any problem and allow a reduction of the sample volume up to 50 %.

- Cell growth area: 15.0 mm²
- Working volume: 15 – 175 μ l



Bulk Packaging

For selected products Greiner Bio-One also offers user-friendly bulk packaging (Fig. 1). Additional products are available in bulk pack upon request.





Figure 1: User-friendly bulk packaging




96 Well Polystyrene Cell Culture / Suspension Culture Microplates

solid bottom, clear


- 96 Well Microplates p. 2 | 7
- Sealers, Lids and CapMats p. 12 | 2 ff.
- Barcode Labelling p. 14 | 4

Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic




Cat.-No.	650 160	650 180	650 185	651 160	651 180
Well profile	U-bottom	U-bottom	U-bottom	V-bottom	V-bottom
Bottom	solid	solid	solid	solid	solid
Colour	clear	clear	clear	clear	clear
Growth area per well [mm²]	35	35	-	28	28
Working volume per well [µl]	40 – 280	40 – 280	40 – 280	40 – 200	40 – 200
TC surface treatment	+	+	-/suspension	+	+
Sterile	+	+	+	+	+
Lid	-	+	+	-	+
Quantity per bag/case	1/100	1/100	1/60	1/100	1/100



Cat.-No.	655 160	655 162	655 180	655 182	655 185
Well profile	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well
Bottom	solid	solid	solid	solid	solid
Colour	clear	clear	clear	clear	clear
Growth area per well [mm²]	34	34	34	34	-
Working volume per well [µl]	25 – 340	25 – 340	25 – 340	25 – 340	25 – 340
TC surface treatment	+	+	+	+	-/suspension
Sterile	+	+	+	+	+
Lid	-	-	+))	+))	+))
Quantity per bag/case	1/100	5/100	1/100	10/160	1/60

+) with condensation rings




96 Well Polystyrene Cell Culture Microplates

solid bottom, white / black

- 96 Well Microplates p. 2 | 8
- Sealers, Lids and CapMats p. 12 | 2 ff.
- Barcode Labelling p. 14 | 4

Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic




Cat.-No.	655 073	655 083	655 079	655 086
Well profile	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well
Bottom	solid	solid	solid	solid
Colour	white	white	black	black
Growth area per well [mm²]	34	34	34	34
Working volume per well [µl]	25 – 340	25 – 340	25 – 340	25 – 340
TC surface treatment	+	+	+	+
Sterile	+	+	+	+
Lid	-	+*)	-	+*)
Quantity per bag/case	10/40	8/32	10/40	8/32

*) with condensation rings

655 087
655 090655 088
655 098**96 Well Polystyrene Cell Culture Microplates****µClear®, white / black**

- 96 Well Microplates p. 2 | 8
- Sealers, Lids and CapMats p. 12 | 2 ff.
- Barcode Labelling p. 14 | 4

Free of detectable
DNase, RNase,
human DNA
non-pyrogenicnon-
cytotoxic

Cat.-No.	655 088	655 098	655 087	655 090
Well profile	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well	F-bottom/ chimney well
Bottom	µClear®	µClear®	µClear®	µClear®
Colour	white	white	black	black
Growth area per well [mm²]	34	34	34	34
Working volume per well [µl]	25 – 340	25 – 340	25 – 340	25 – 340
TC surface treatment / Sterile	+ / +	+ / +	+ / +	+ / +
Lid	-	+ ^{*)}	-	+ ^{*)}
Quantity per bag/case	10/40	8/32	10/40	8/32

*) with condensation rings

2 HTS-
Microplates3 Immunology/
HLA4 Microbiology/
Bacteriology5 Tubes/Multi-
Purpose Beakers6 Liquid
Handling7 Molecular
Biology8 Protein
Crystallisation

9 Separation

10 Biochips/
Microfluidics11 Cryo-
Techniques12 Lids/Sealers/
CapMats13 Reaction Tubes/
Analyser Cups


14 Accessories

675 090

675 098

96 Well Half Area Polystyrene Cell Culture Microplates**solid bottom, clear / white / black****µClear®, white / black**

- 96 Well Half Area Microplates p. 2 | 9

Free of detectable
DNase, RNase,
human DNA
non-pyrogenicnon-
cytotoxic

Cat.-No.	675 180	675 083	675 086	675 098	675 090
Well profile	half area	half area	half area	half area	half area
Bottom	solid	solid	solid	µClear®	µClear®
Colour	clear	white	black	white	black
Growth area per well [mm²]	15	15	15	15	15
Working volume per well [µl]	15 – 175	15 – 175	15 – 175	15 – 175	15 – 175
TC surface treatment / Sterile	+ / +	+ / +	+ / +	+ / +	+ / +
Lid	+	+	+	+	+
Quantity per bag/case	8/32	8/32	8/32	8/32	8/32