
Microscale makes "macro-sense"



· Reduced chemical and disposal costs


· Use with Williamson textbook: "Macroscale and Microscale Organic Experiments"




Kits are packed in a rugged polyethylene storage case (2 1/2" H x 13 1/2" W x 10" D) with a die cut, foam insert. They fit easily into lab bench drawers.

Patented Thermoplastic Connectors


Fast, easy-to-use, positive connections without threads or standard taper joints. US Patent #4,803,053





For additional kits, see our website: www.kimble-chase.com

Access



Exclusive to Kimble Chase!

Williamson Microscale Kits

All the components necessary to perform the Williamson experiments described in Macroscale and Microscale Organic Experiments are included.

Williamson Microscale Organic Chemistry Kit

Improved Standard Microscale Kit with Thermoplastic Connectors

These MICROFLEX® Kits feature patented chaste connectors and specially tooled glassware components providing great versatility. If replacement parts are needed, individual components may be ordered separately.

Catalog No.	Description		2012 Pric
748100-0000	Williamson Improved Standard Microscale	e Kit	266.12
			•
Kit Compo	onents		1
Catalog No.	Description	Item	Kit Qty.
748001-1000	Connecting Adapter	Α	1
748002-0000	Filter Adapter, Pluro Stopper	В	1
748003-1000	Thermometer Adapter	С	1
411800-1015	Centrifuge Tube, 15 mL	D	1
420166-1001	Chromatography Column,	Ε	1
	Glass 150 mm w/Disc		
420163-1500	One Way Stopcock, PE/Nylon	F	1
748005-1000	Connector w/Support Rod	G	1
748005-1001	Connector Only	H	1
748006-1100	Distillation Column	l.	1
748007-1060	Distillation Head, 60 mm	J	1
748008-1005	Short Neck Flask, 5 mL	K	1
748009-1005	Long Neck Flask, 5 mL	L	1
882300-0010	Erlenmeyer Flask, 10 mL	M	3
748010-1025	Filter Flask, 25 mL	N	1
748012-1000	Hirsch Funnel, with Disc	0	1
748013-1000	Reaction Tube, 10 x 100 mm	Р	5
748014-1000	Spatula, S.S.	Q	1
774261-0008	Sleeve Stoppers	R	2
748019-0001	Syringe, Polyethylene	S	1
420823-0116	Tubing, PTFE, 1/16" for 20 gauge needle	T	2
764100-1001	Pipet, 1 mL	U	1

Stir Bars, 4x12 mm, 2 per vial

Connector Only, Viton®

791145-0021


748005-1651

Limited Time Offer! Purchase Five 748100-0000 Kits and Get an Additional Kit of the Same Type for FREE. See details on back.



888-546-2531 info@kimble-chase.com www.kimble-chase.com

Williamson Standard Kit, 748000-0000



F -2.5 -1.5 -1.5 -1.5 -1.5 -1.5 -1.5 -1.5 -1
13 15 15 15 15 15 15 15 15 15 15 15 15 15
S P
G J K
T

Williamson Intermediate Kit, 748110-0000

Catalog No.	Description	2012 Price
748000-0000	Williamson Standard Microscale Kit	257.11


(it Components			

Kit Compo	nents		
Catalog No.	Description	Item	Kit Qty
748001-1000	Connecting Adapter	Α	1
748002-0000	Filter Adapter, Pluro Stopper	В	1
748003-1000	Thermometer Adapter	С	1
411800-1015	Centrifuge Tube, 15 mL	D	1
420166-1001	Chromatography Column, Glass 150 mm w/Disc	Ε	1
420163-1500	One Way Stopcock, PE/Nylon	F	1
748005-1651	Connector Only, Viton®	G	1
748005-1000	Connector w/Support Rod	Н	1
748006-1100	Distillation Column	I	1
748007-1060	Distillation Head, 60 mm	J	1
748008-1005	Short Neck Flask, 5 mL	K	1
748009-1005	Long Neck Flask, 5 mL	L	1
882300-0010	Erlenmeyer Flask, 10 mL	M	4
748010-1025	Filter Flask, 25 mL	N	1
748012-1000	Hirsch Funnel, with Disc	0	1
748013-1000	Reaction Tube, 10 x 100 mm	Р	4
748014-1000	Spatula, S.S.	Q	1
774261-0008	Sleeve Stoppers	R	2
748019-0001	Syringe, Polyethylene	S T	1
420823-0116	Tubing, PTFE,1/16" fits 20 gauge needle	T	2

Catalog No.	Description	2012 Price
748110-0000	Williamson Intermediate Microscale Kit	134.95

Kit Components

Catalog No.	Description	Item	Kit Qty
748001-1000	Connecting Adapter	Α	1
748002-0000	Filter Adapter, Pluro Stopper	В	1
748003-1000	Thermometer Adapter	С	1
411800-1015	Centrifuge Tube, 15 mL	D	1
420166-1001	Chromatography Column, Glass 150 mm w/Disc	Ε	1
420163-1500	One Way Stopcock, PE/Nylon	F	1
748005-1000	Connector w/Support Rod	G	1
748005-1651	Connector Only, 'Viton®	Н	1
748006-1100	Distillation Column		1
748007-1060	Distillation Head, 60 mm	J	1
748008-1005	Short Neck Flask, 5 mL	K	1
748010-1025	Filter Flask, 25 mL	N	1
748012-1000	Hirsch Funnel, with Disc	0	1
748013-1000	Reaction Tube, 10 x 100 mm	Р	2
748019-0001	Syringe, Polyethylene	S	1
420823-0116	Tubing, PTFE, 1/16" fits 20 gauge needle	T	2







Williamson Macroscale Organic Chemistry Kit


KEM-KIT® with \$ 14/20 Joints

Designed to address the experiments detailed in the Williamson text Macroscale and Microscale Organic Experiments. These experiments can also be performed with the 748000 Series Microscale Kits.

Students use the separatory funnel for extractions and for liquid addition in reaction assemblies. Classical atmospheric, steam, fractional, and vacuum distillations are possible with kit components. Apparatus set-up time is minimal, providing the student with more lab time. Smaller quantities afford a greater degree of safety, improve air quality in the lab, and substantially reduced disposal costs.

Catalog No.	Description		2012 Price	
269330-0000	Williamson Organic Chemistry Macro	scale Kit	466.29	
Kit Components				
Catalog No.	Description	Item	Kit Qty.	
251450-0015	Conical Flask, 15 mL	A	1	

Catalog No.	Description	Item	Kit Qty.
251450-0015	Conical Flask, 15 mL	А	1
273755-0000	Claisen Adapter	В	1
274960-0000	Connecting Adapter	С	1
276405-0000	Vacuum Adapter	D	1
282555-0000	Water Jacketed Condenser	E	1
286100-0000	Distilling Column	F	1
288260-0001	Distillation Receiver	G	1
294300-0050	Kjeldahl Shaped Flask, 50 mL	Н	1
294300-0100	Kjeldahl Shaped Flask, 100 mL	I	1
298281-0125	Separatory Funnel, 125 mL	J	1
411800-0015	Centrifuge Tube	K	1
420163-1500	Stopcock, PE/Nylon	L	1
420166-1001	Chromatography Column	M	1
675300-0014	Polyacetal Clamp, Size 14	N	4
748002-0000	Filter Adapter, Pluro Stopper	0	1
748010-1025	Filter Flask, 25 mL	Р	1
748012-1000	Hirsch Funnel	Q	1
850400-1420	Hex Head Stopper, \$ 14/20	R	1





To redeem this offer, please fax a copy of your completed coupon along with proof of purchase (invoice or P.O. dated between 1/1/12 and 12/31/12) to FREE Kimble Chase Microscale Kit at 865-354-3823.

Name: Company:

Address:

City: State / Province:

Postal Code: Country:

Email: Phone: Fax:

Redemption must be received by 1/15/13. Limited to two redemptions per customer. Please allow 6-8 weeks for delivery. International participant redemptions to be fulfilled by their local office and subject to product substitution and extension of delivery time. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. Kimble Chase distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Trademarks contained herein are the property of their respective owners.

Printed in USA. Trademarks are the property of their respective owners. Specifications and prices subject to change without notice. Prices shown are US dollars. PDF0412KC-1105R

See our complete line of laboratory glassware on the web:

www.kimble-chase.com info@kimble-chase.com

Your Source for:

Kimble, Kontes and Chase Brands

Call: 888-546-2531 Fax: 856-794-9762



