

KIMAX-HQ™ NMR Tubes processed with our proprietary RAY-SORB® treatment.

RAY-SORB® protects your valuable samples from the harmful effect of visible and UV wavelengths, while keeping all the critical characteristics of your specific tube needs intact.

RAY-SORB® is also available as a custom option for many other NMR items.


Call or visit our website for more information on our new line of RAY-SORB® NMR tubes from Kimble Chase®.

Compatible with instruments from the following manufacturers.

	APPEN			
•	R	ru	ke	r
			$\sim$	

- Agilent
- · Jeol and many others.

	h			
		5mr		
<b>Catalog No.</b> 897241-9001	OD	Length		Case (
897241-9001	5mm	8"	600	5
897245-9001	5mm	8"	800	5



For additional product, technical, or distributor information, please visit www.kimble-chase.com or call 800.451.4351. Customers outside the United States, please call: +1.865.717.2600 Ext. 1 or contact your local Kimble Chase® sales office listed below.

NORTH AMERICA Tel: 800.451.4351 Fax: 865.354.3853

EUROPE Tel: 49.3693.8136.0 Fax: 49.3693.8136.33 LATIN AMERICA Tel: 52.442.290.7446 Fax: 52.442.290.7442

ASIA / PACIFIC Tel: 86.10.8739.5506, .5509, .1956 Fax: 86.10.8739.1990

KIMBLE CHASE

234 Cardiff Valley Road, Rockwood, TN 37854 USA www.kimble-chase.com