

Your **Power** for Health



greiner bio-one



**VACUETTE®**


**SAFETY Blood  
Collection Sets**

Our Innovations  
for Your Safety

[www.gbo.com/preanalytics](http://www.gbo.com/preanalytics)

# SAFETY Blood Collection Sets

Specially developed for blood collection from patients with difficult vein conditions. Protect yourself and your employees from the risk of needlestick injury. Every needlestick injury is one too many.



## Activation

Complete activation of the safety mechanism is clearly indicated by an audible "click".

## Safety Mechanism

Simple activation of the safety mechanism via lightly pressing in both sides of the hub.

## Transparent View Window

The transparent view window provides clear vein entry indication.

## Especially Safe

Particularly safe due to activation of the safety mechanism while the needle is still in the vein.

## More Flexibility

Assorted tube lengths offer more flexibility for blood collection.

## Optional

use as an infusion set.

Item No. 450090-92

# SAFETY Blood Collection Sets Combination Products

As a practical combination product, the SAFETY Blood Collection Set is also available pre-assembled with:

## SAFETY Blood Collection Set + Luer Adapter



Item No. 450081-84,  
Item No. 450095-96,  
Item No. 450099

- Pre-attached set with Luer adapter in sterile blister pack.
- Handling is optimised.

## SAFETY Blood Collection Set + Holder



Item No. 450085-87,  
Item No. 450080

- Pre-attached set with holder in sterile blister pack.
- No threading necessary.
- A practical, ready-to-use solution that saves time.


## SAFETY Blood Collection Set + Blood Culture Holder



Item No. 450182-85

- Pre-attached set with Blood Culture Holder in sterile blister pack.
- Designed for use with most commonly used blood culture bottles and **VACUETTE**<sup>®</sup> blood collection tubes in one procedure.
- The set is sterile and pre-assembled for convenience, cost effectiveness and ensures your patient only has to undergo one venipuncture for all procedures.

## Instructions for Use \*




**1** Remove protective cap




**2** Venipuncture




**3** Insert the tube




**4** Remove the last tube



**5** Activate the safety mechanism (by pressing down and pulling back the lighter coloured area)



**6** Audible click indicates activation



**7** Activated SAFETY Blood Collection Set



**8** Disposal

\* For the complete handling recommendation please refer to the Instructions for Use for each product.

## SAFETY Blood Collection Sets

Item. No.	Description	Size	Packaging		
			Inner	Outer	
<b>SAFETY Blood Collection/Infusion Sets (sterile)</b>					
450090	SAFETY Blood Collection/Infusion Set	tubing length: 19cm	25G x 3/4"	50 pcs.	1.000 pcs.
450091	SAFETY Blood Collection/Infusion Set	tubing length: 19cm	21G x 3/4"	50 pcs.	1.000 pcs.
450092	SAFETY Blood Collection/Infusion Set	tubing length: 19cm	23G x 3/4"	50 pcs.	1.000 pcs.
<b>SAFETY Blood Collection Sets + Luer Adapters (sterile)</b>					
450083	SAFETY Blood Collection Set + Luer Adapter	tubing length: 10cm	21G x 3/4"	50 pcs.	1.000 pcs.
450081	SAFETY Blood Collection Set + Luer Adapter	tubing length: 19cm	21G x 3/4"	50 pcs.	1.000 pcs.
450095	SAFETY Blood Collection Set + Luer Adapter	tubing length: 30cm	21G x 3/4"	50 pcs.	1.000 pcs.
450084	SAFETY Blood Collection Set + Luer Adapter	tubing length: 10cm	23G x 3/4"	50 pcs.	1.000 pcs.
450082	SAFETY Blood Collection Set + Luer Adapter	tubing length: 19cm	23G x 3/4"	50 pcs.	1.000 pcs.
450096	SAFETY Blood Collection Set + Luer Adapter	tubing length: 30cm	23G x 3/4"	50 pcs.	1.000 pcs.
450099	SAFETY Blood Collection Set + Luer Adapter	tubing length: 30cm	25G x 3/4"	50 pcs.	1.000 pcs.
<b>SAFETY Blood Collection Sets + Holders (sterile)</b>					
450087	SAFETY Blood Collection Set + Holder	tubing length: 10cm	21G x 3/4"	24 pcs.	240 pcs.
450085	SAFETY Blood Collection Set + Holder	tubing length: 19cm	21G x 3/4"	24 pcs.	240 pcs.
450086	SAFETY Blood Collection Set + Holder	tubing length: 19cm	23G x 3/4"	24 pcs.	240 pcs.
450080	SAFETY Blood Collection Set + Holder	tubing length: 30cm	25G x 3/4"	24 pcs.	240 pcs.
<b>SAFETY Blood Collection Sets + Blood Culture Holders (sterile)</b>					
450182	SAFETY Blood Collection Set + Blood Culture Holder	tubing length: 19cm	21G x 3/4"	24 pcs.	240 pcs.
450183	SAFETY Blood Collection Set + Blood Culture Holder	tubing length: 19cm	23G x 3/4"	24 pcs.	240 pcs.
450184	SAFETY Blood Collection Set + Blood Culture Holder	tubing length: 30cm	21G x 3/4"	24 pcs.	240 pcs.
450185	SAFETY Blood Collection Set + Blood Culture Holder	tubing length: 30cm	23G x 3/4"	24 pcs.	240 pcs.

These devices have no components made of natural dry rubber.



For further information, please visit our website [www.gbo.com/preanalytics](http://www.gbo.com/preanalytics) or contact us:

**Austria (Main Office)**

Greiner Bio-One GmbH  
Phone +43 7583 6791-0  
Fax +43 7583 6318  
E-Mail [office@at.gbo.com](mailto:office@at.gbo.com)

**Germany**

Greiner Bio-One GmbH/Preanalytics  
Tel +49 7022 948-0  
Fax +49 7022 948-514  
E-Mail [office@de.gbo.com](mailto:office@de.gbo.com)

**Spain**

VACUETTE Espana S.A.  
Phone +34 91 652 77 07  
Fax +34 91 652 33 35  
E-Mail [info@vacuette.es](mailto:info@vacuette.es)

**United Kingdom**

Greiner Bio-One Ltd.  
Phone +44 1453 8252 55  
Fax +44 1453 8262 66  
E-Mail [info@uk.gbo.com](mailto:info@uk.gbo.com)

**Brazil**

Greiner Bio-One Brasil  
Phone +55 19 3468-9600  
Fax +55 19 3468-9621  
E-Mail [office@br.gbo.com](mailto:office@br.gbo.com)

**Hungary**

Greiner Bio-One Hungary Kft.  
Phone +36 96 213 088  
Fax +36 96 213 198  
E-Mail [office@hu.gbo.com](mailto:office@hu.gbo.com)

**Switzerland**

Greiner Bio-One VACUETTE  
Schweiz GmbH  
Phone +41 7 12 28 55 22  
Fax +41 7 12 28 55 21  
E-Mail [office@ch.gbo.com](mailto:office@ch.gbo.com)

**USA**

Greiner Bio-One North America Inc.  
Phone +1 704 261-7800  
Fax +1 704 261-7899  
E-Mail [info@us.gbo.com](mailto:info@us.gbo.com)

**China**

Greiner Bio-One Suns Co., Ltd.  
Phone +86 10 83 55 19 91  
Fax +86 10 63 56 69 00  
E-Mail [office@cn.gbo.com](mailto:office@cn.gbo.com)

**India**

Greiner Bio-One INDIA Pvt., Ltd.  
Phone +91 120 456 8787  
Fax +91 120 456 8788  
E-Mail [info@gboindia.com](mailto:info@gboindia.com)

**Thailand**

Greiner Bio-One Thailand Ltd  
Phone +66 38 4656 33  
Fax +66 38 4656 36  
E-Mail [office@th.gbo.com](mailto:office@th.gbo.com)

**France**

Greiner Bio-One SAS  
Phone +33 1 69 86 25 25  
Fax +33 1 69 86 25 35  
E-Mail [office@fr.gbo.com](mailto:office@fr.gbo.com)

**Netherlands**

Greiner Bio-One B.V.  
Phone +31 172 4209 00  
Fax +31 172 4438 01  
E-Mail [info@nl.gbo.com](mailto:info@nl.gbo.com)