


Locking Dialysis Membrane Clamps


Superior in sealing dialysis tubing, our clamps stay locked until you release them, and they do not float.

Constructed of Polyamide 6.6 from DuPont (NYLON). Available in three sizes to accommodate flat widths up to 46 mm, 65 mm or 105 mm. Ten clamps per package.



Part No.	Max Sealing Width	Total Length
CB-10100	105mm	130mm
CB-1070	65mm	90mm
CB-1050	46mm	65mm

FILLER FUNNELS




Cellu•Sep Filler Funnels make loading and unloading dialysis tubing easy: Cut membrane tube to length and clamp at one end, slip the open end onto the polypropylene funnel and secure with the rubber loop fastener. Open snap cap, fill with sample, and hang in container. Funnels are supported by a three-notched universal hanger that fits over the edge of common vessels.

MACRO FILLER FUNNEL *Part No. MFPL-05*

For 23 mm, 25 mm, 32 mm, 40 mm and 45 mm flat widths. Package of five funnels with rubber loops.

MINI FILLER FUNNEL *Part No. MFPS-05*

For 10 mm and 19 mm flat widths. Package of five funnels and rubber loops.



BAG HOLDERS

Part No. MBHS-05

These strong polypropylene holders suspend dialysis bags in containers. Easy to mark for identifying sequence when dialyzing multiple samples. Package of 5 holders.