



## Heavy Duty Polypropylene Vacuum Bottles

Roush Life Sciences Heavy Duty Vacuum Bottles combine maximum ease of use and versatility, with the highest level of strength and durability.

- ❖ Heavy Duty Vacuum bottles have thick walls that can hold full vacuum: 14.8PSI/765mmHg for 24 hours, making them excellent replacements for glass HPLC containers, and ideal for applications involving liquid aspiration, cell harvesting, and degasification.
- ❖ Heavy Duty Vacuum Bottles are compatible with 50mm (1 and 2 L bottles) and 80mm (4L) VersaCaps™ and adapters, enabling easy to customize the bottle to fit your application.
- ❖ All Bottles includes Cap with closed and 1/4" HB adapters
- ❖ Bottle material and volume are molded into the side of the bottle for easy identification, along with metric graduation marks certified accurate to ±5%.
- ❖ Easy opening & closing without the need to disconnect tubes
- ❖ Polypropylene Heavy Duty Vacuum Bottles, as well as all caps and adapters are autoclavable, at 15 PSI and 121°C for 15 minute s

Cat. #	Bottles Size	dimensions
165-1102	1 liter	11x11x28 cm
165-2102	2 liters	13x13x33 cm
165-4202	4 liters	18x18x38 cm



Other materials (PC, PETG, HDPE) upon request



August 2012



לפרטים נוספים:  
 דורית: 054-2209481  
 אינה: 054-4949145  
 זהר: 054-6541442  
 עדי: 054-6426636  
 דב: 054-6455268