

Cryo.s™ Rack Scanner

Ideal for Labs Scanning Frozen Racks

With the Cryo.s™ Rack Scanner Greiner Bio-One provides an optimum solution for scanning ambient temperature or frozen Cryo.s[™] with Datamatrix and Cryo.s[™] Biobanking tubes in SBS-formatted racks.

- Scans and decodes one rack in 8-10 seconds 6
- Insulates rack from scanner surface 6
- Anti-fog treatment on inside scanning surface 6
- Internal temperature control system can be turned on or off as needed 6
- Single tube scan feature 6
- Decodes all SBS rack formats 6
- Easy one-touch button for initiation of scan 6
- 2D embedded barcode for Rack ID capture 6
- Small footprint 6
- Integration friendly design 6
- Auto-rescan feature improves decoding results 6
- 6 5 year warranty

The scanner comes with a comprehensive software package comprising the following functions:

- Scanning of 2D barcodes on tubes and racks, output of decoded barcode and rack ID 6
- Rescan or manually edit barcode data 6
- Single tube scan displays 2D code on screen and pastes into any open application 6
- Project based workflow to allow easy switching between rack formats 6
- Data output files are customizable .csv format 6
- SDK and ActiveX toolkits are included to enable transfer of output files directly to sample 6 management software, LIMS or to integrate with other laboratory robotics
- Lifetime upgrades guarantee for software 6

The scanner comes with a comprehensive software package comprising the following functions:


- PC; Windows 10, 8, 7, Vista 6
- Two USB ports 6

Germany (Main office): Greiner Bio-One GmbH, info@de.gbo.com | Austria: Greiner Bio-One GmbH, office@at.gbo.com Belgium: Greiner Bio-One BVBA/SPRL, info@be.gbo.com | Brazil: Greiner Bio-One Brasil, office@br.gbo.com | China: Greiner Bio-One Suns Co. Ltd., office@cn.gbo.com France: Greiner Bio-One SAS, infos@fr.gbo.com | Japan: Greiner Bio-One Co. Ltd., info@jp.gbo.com | Netherlands: Greiner Bio-One B.V., info@nl.gbo.com UK: Greiner Bio-One Ltd., info@uk.gbo.com | USA: Greiner Bio-One North America Inc., info@us.gbo.com WW



Cryo.s[™] Rack Scanner

Ideal for Labs Scanning Frozen Racks



Ordering Information		
Cat. No.	Product Description	Box contents
849 050	Cryo.s™ Rack Scanner, US edition	1 scanner, power cables, USB cables, drivers and manual
849 070	Cryo.s™ Rack Scanner, EU/UK/Asia edition	1 scanner, power cables, USB cables, drivers and manual

*A sample management software is not part of the included software package.

Germany (Main office): Greiner Bio-One GmbH, info@de.gbo.com | Austria: Greiner Bio-One GmbH, office@at.gbo.com Belgium: Greiner Bio-One BVBA/SPRL, info@be.gbo.com | Brazil: Greiner Bio-One Brasil, office@br.gbo.com | China: Greiner Bio-One Suns Co. Ltd., office@cn.gbo.com France: Greiner Bio-One SAS, infos@fr.gbo.com | Japan: Greiner Bio-One Co. Ltd., info@jp.gbo.com | Netherlands: Greiner Bio-One B.V., info@nl.gbo.com UK: Greiner Bio-One Ltd., info@uk.gbo.com | USA: Greiner Bio-One North America Inc., info@us.gbo.com WWW