




SealPlate® RollMate™



[Click to download artwork](#)

www.ExcelScientific.com/DS/

SealPlate® RollMate™: pre-scored film in a convenient, easy to use roll format

Application:

Ideal for ELISA, incubation and storage

Temperature Range:

-40°C to +120°C

Catalog #	Description	Quantity
SP2-RL-500	SealPlate RollMate	500/Roll

Features and Benefits

- ▶ An economical, efficient way to rapidly seal microplates; faster than applying single sheeted film
- ▶ Self-standing, non-slip supports hold pre-scored film roll securely in place for quick and easy dispensing
- ▶ Continuous backing liner can be torn off or fed directly into a waste basket for disposal
- ▶ Roll tension can be easily adjusted to user preference (e.g. loose, tight)
- ▶ SealPlate film minimizes evaporation, reduces cross-contamination and prevents spillage
- ▶ Eliminates the “edge effect”, excess evaporation in the outside wells, which can occur with other films and plate lids
- ▶ Useful as a simple cover to protect samples from external contamination, dust and debris
- ▶ 2.0mil rigid polyester film features good clarity and easy handling
- ▶ Single, writable end-tab for easy application and plate identification
- ▶ High integrity packaging prevents damage during transit; poly bag eliminates dust and particulate contamination

Bundling Opportunities

- ▶ ELISA plates
- ▶ ELISA reagents and buffers
- ▶ Plate washers
- ▶ Plate readers
- ▶ Pipettors
- ▶ Pipet tips